


AutoCart™

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Reference for AutoCart

1. Albrecht F. et al. Closure of Osteochondral Lesions Using Chondral Fragments and Fibrin Adhesive. Arch Orthop Trauma Surg (1983) 101: 213 - 217
2. Stone K. Articular Cartilage Paste Grafting to Full-Thickness Articular Cartilage Knee Joint Lesions: A 2- to 12-Year Follow-Up, Arthroscopy: The Journal of Arthroscopic and Related Surgery, Vol 22, No 3 (March), 2006: pp 291 – 299
3. Christensen et al. Autologous Dual-Tissue Transplantation for Osteochondral Repair: Early Clinical and Radiological Results; Cartilage 2015, Vol. 6 (3) 166 - 173
4. Massen F. et al. One-Step Autologous Minced Cartilage Procedure for the Treatment of Knee Joint Chondral and Osteochondral Lesions; The Orthopaedic Journal of Sports Medicine, 7(6),
5. Feeney et al. Autologous Cartilage Particulate for Treatment of Cartilage Defects: Impact of Different Arthroscopic Shavers on Viability and In Vitro Migration; ORS poster 2020
6. Olena Virchenko et al: Independent and additive stimulation of tendon repair by thrombin and platelets, Acta Orthopaedica 2006; 77 (6): 960–966.
7. Kosuke Uehara et al: Effect of Fibrin Formulation on Initial Strength of Tendon Repair and Migration of Bone Marrow Stromal Cells in Vitro. J Bone Joint Surg Am. 2015;97:1792-8

Reference for Shaver Blades

1. Lu Y. et al. Minced Cartilage without Cell Culture Serves as an Effective Intraoperative Cell Source for Cartilage Repair. Journal of Orthopaedic Research June 2006: 1261 - 1270
2. Feeney et al. Autologous Cartilage Particulate for Treatment of Cartilage Defects: Impact of Different Arthroscopic Shavers on Viability and In Vitro Migration; ORS poster 2020

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