IOBP[®]

IntraOsseous BioPlasty[®]

Arthrex[®] 

Table of Contents

IntraOsseous BioPlasty®	03
IOBP® Decompression Kits	04
Angel™ Bone Marrow Concentration	05
Thrombinator™ System.....	06
Preparation of Injection Solution	07
自家トロンビン液の調整方法.....	08
BMCと自家トロンビン液の混和.....	10
Surgical Technique	11
IntraOsseous BioPlasty ガイドを使用する場合	12
IntraOsseous BioPlasty デコンプレッションデバイスを使用する場合	13
Ordering Information.....	15
References	16

IntraOsseous BioPlasty®

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Reference for Angel Bone Marrow Concentration

1. Mark R Nazal et al: Connective Tissue Progenitor Analysis of Bone Marrow Aspirate Concentrate Harvested From the Body of the Ilium During Arthroscopic Acetabular Labral Repair. Arthroscopy, 2020 May;36(5):1311-1320.
2. Hisashi Sugaya et al: Comparative Analysis of Cellular and Growth Factor Composition in Bone Marrow Aspirate Concentrate and Platelet-Rich Plasma. Bone Marrow Res, 2018 Feb 25;2018:1549826
3. Christof Pabinger et al: Intra-articular injection of bone marrow aspirate concentrate (mesenchymal stem cells) in KL grade III and IV knee osteoarthritis: 4 year results of 37 knees. Scientific Reports, (2024) 14:2665
4. Christopher Centeno et al: A Randomized Controlled Trial of the Treatment of Rotator Cuff Tears with Bone Marrow Concentrate and Platelet Products Compared to Exercise Therapy: A Midterm Analysis. Stem Cells Int. 2020 Jan 30;2020:5962354
5. Francesca Vannini et al: Bone marrow aspirate concentrate and scaffold for osteochondral lesions of the talus in ankle osteoarthritis: satisfactory clinical outcome at 10 years. Knee Surgery, Sports Traumatology, Arthroscopy (2021) 29:2504–2510
6. Philippe Hernigou et al: Subchondral bone or intra-articular injection of bone marrow concentrate mesenchymal stem cells in bilateral knee osteoarthritis: what better postpone knee arthroplasty at fifteen years? A randomized study, Int Orthop, 2021 Feb;45(2):391-399.
7. Philippe Hernigou et al: Human bone marrow mesenchymal stem cell injection in subchondral lesions of knee osteoarthritis: a prospective randomized study versus contralateral arthroplasty at a mean fifteen year follow-up. Int Orthop. 2021 Feb;45(2):365-373.

Reference for Thrombinator System

1. Olena Virchenko et al: Independent and additive stimulation of tendon repair by thrombin and platelets, Acta Orthopaedica 2006; 77 (6): 960–966.
2. Kosuke Uehara et al: Effect of Fibrin Formulation on Initial Strength of Tendon Repair and Migration of Bone Marrow Stromal Cells in Vitro. J Bone Joint Surg Am. 2015;97:1792-8

Reference for IntraOsseous BioPlasty

1. Fuchs A et al.: Implant Augmentation for Trochanteric Fractures with an Innovative, Ready to Use Calcium-Phosphate-Cement. Journal of Orthopaedics and Bone Research. 2019
2. Niall DM, Bobic VB. Bone Bruising and Bone Marrow Edema Syndromes: Incidental Radiological Findings or Harbingers of Future Joint Degeneration? International Society of Arthroscopy, Knee Surgery, and Orthopaedic Sports Medicine. Retrieved from: <https://www.isakos.com/innovations/niall>
3. Smith RC, Constable RT, Reinhold C, McCauley T, Lange RC, McCarthy S. Fast Spin Echo STIR Imaging. Journal of Computer Assisted Tomography. 1994; 18(2):209-213

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